


INOVO SENSE






*Sense Of
Confidence*

THE INOVOSENSE® FAMILY

	SHINE	SILK	SMILE	STYLE	SHAPE
Product					
Cross-Linking Level	Non cross-linked	●○○○	●●○○	●●●○	●●●●
Gel Texture					
G'	70	200	300	500	550
HA Concentration	10 mg / ml	19 mg / ml	20 mg / ml	23 mg / ml	25 mg / ml
Lidocaine	0.3 %	0.3 %	0.3 %	0.3 %	0.3 %
Needle Size	30 G ¹ / ₂	30 G ¹ / ₂	30 G ¹ / ₂ 27 G ¹ / ₂	27 G ¹ / ₂	27 G ¹ / ₂ 25 G ¹ / ₂
Volume	1 ml	2 × 1 ml	2 × 1 ml	2 × 1 ml	2 × 1 ml
Indications	 <ul style="list-style-type: none"> ● Skin hydration for: <ul style="list-style-type: none"> -Face -Neck -Décolletage -Back of hands 	 <ul style="list-style-type: none"> ● Correction of superficial wrinkles, such as: <ul style="list-style-type: none"> -Crow's feet -Glabellar lines -Perioral lines ● Lip contouring 	 <ul style="list-style-type: none"> ● Correction of moderate wrinkles, such as: <ul style="list-style-type: none"> -Nasolabial wrinkles -Marionette lines -Oral commissures ● Lip volume 	 <ul style="list-style-type: none"> ● Correction of severe wrinkles, such as: <ul style="list-style-type: none"> -Nasolabial folds ● Facial lifting in: <ul style="list-style-type: none"> -Cheeks and cheekbones -Chin -Temples 	 <ul style="list-style-type: none"> ● Shape and restoration of facial volume in: <ul style="list-style-type: none"> -Cheeks and cheekbones -Chin -Jawline -Temples
Depth of Injection	Epidermis or dermo-epidermal junction	Superficial to mid Dermis	Mid to deep Dermis	Deep dermis or subcutaneous	Deep dermis or subcutaneous



 +98 21 82430 - 752
 info@medisaara.com
 inovosense.com
 inovosense
 No.21, Jooybar Alley, Shad St., Mollasadra Ave., Vanak Sq., Tehran, Iran

INOVO
SENSE