



■ Fiorage Family



□ FIORAGE S


Dosage:

16 mg Cross-linked hyaluronic acid
Each syringe contains 1 mL of **Fiorage S**

Indication:

injection into the middle layer of the skin
fill wrinkles and superficial pits
shaping the lips






Dosage:

20 mg Cross-linked Hyaluronic acid gel
Each syringe contains 1 ml of **Fiorage M**

Indications:

natural correction of superficial lines and wrinkles
the deepest of wrinkles away from the mouth
Lip volume augmentation
restoration of volume in the back of the hand





□ FIORAGE L


Dosage:

24 mg Cross-linked hyaluronic acid
Each syringe contains 1 mL of **Fiorage L**

Indications:

injection into the middle layer of the skin
fill wrinkles and pits, restoration of volume
in the back of the hand





Dosage:

26 mg Cross-linked hyaluronic acid
Each syringe contains 1 mL of **Fiorage XL**

Indications:

injection into the deeper layers of the skin
repair or enhance the fullness of the face
(cheek, chin, jawline sculpting)

